



# OPERATIONSLAMPA POWER LED II

## EN BEKVÄM PRISVÄRD OPERATIONSLAMPA

Power LED II kombinerar den indirekta teknologin som konsolideras av PentaLED-seriens framgång med Unicas eleganta och kompakta design.

POWER LED II är en högpresterande lampa, prisvärd i sin klass och smart konstruerad, för alla olika kirurgers behov.



## SPECIFIKATIONER

## POWER LED II

Light intensity at 1 m distance (Ec)	160.000 lux
Light head diameter	52 cm
Color temperature	2 selections: 4.500 - 5.000K
Color rendering index (CRI)	96 Ra
Diameter adjustment	Electronic
d10 light field diameter where illuminance reached 10% of Ec	200 mm
Light field diameter adjustable from-to	140 - 280 mm
Depth of illumination IEC 60601-2-41 (L1+L2) at 60%	500 mm
Depth of illumination IEC 60601-2-41 (L1+L2) at 20%	1350 mm
Total radiated energy Ee where the illuminance reaches max level	580 W/m <sup>2</sup>
Ratio between radiated energy Ee and illuminance Ec	3,68
Average led life	> 60.000
Control of the illuminance	20 - 100 %
Electrical absorption	50 W - 60 VA

Direktiv 93/42/EEC - Norm IEC 60601-2-41

### FUNKTIONER



### TILLBEHÖR

